

## INTRODUCTION


B- and T-cell acute lymphocytic leukemia (ALL) affects both children and adults. Cytogenetic and recent molecular genetic findings at diagnosis constitute important, independent prognostic factors in all age groups. Due to the typically poor morphology of the chromosome preparations, conventional cytogenetic analysis can miss one or more recurrent abnormalities and most ALL FISH panels only target select abnormalities. Therefore, a simple, fast, and comprehensive improvement to karyotype analysis would greatly enhance the diagnostic capabilities of clinical service as well as research laboratories for ALL.

## MATERIALS & METHODS

We recently developed and validated a novel molecular technology, Interphase Chromosome Profiling (ICP) (Babu et al., 2017). We tested a variation (ALL-ICP) of our original method for interrogating interphase cells and metaphase chromosomes in ALL samples. The design consisted of a multiplex approach with analysis of six chromosomes per hybridization site. Telomeres and pericentromeric regions on each chromosome are targeted and the resulting fluorescent signals are spectrally distinct and easily recognizable from each other using two dual filter sets from Chroma. An additional probe set targeting common deletions and duplication/amplification was included. Overnight hybridization is done on one slide in four areas and whenever possible 20 metaphase spreads are analyzed for each chromosome. When 20 metaphases are unavailable, interphase nuclei were used to complete the 20-cell analysis. Modified culture conditions included increased time and exposure to dilute concentrations of Colcemid and/or exposure to a chemical known to induce premature chromosome condensation during the last three hours of culture. Forty-eight ALL samples were studied blindly by the laboratory that developed the technology. Three institutions provided cell pellets with known cytogenetic and FISH results.

## RESULTS


Sample Number	Cytogenetic Findings	ICP Results	Additional changes found and/or characterized by ALL-ICP
1	Hyperdiploid, +der	Concordant +	Characterized the derivative
2	Add(5q), inv(7)	Inv missed	Add origin 15q
3	Tetraploid, add(3p), 6(p) 1-3 mar	Concordant +	Deletion 14q, markers origin 13q and 10q
4	Complex abnormalities	Concordant +	14q deletion and dim 14q
5	Complex, add and interstitial deletions	Missed interstitial deletions	Add origin 19q
6	Complex, add(5q), add(17q), interstitial deletions	Concordant +, Missed interstitial deletions	Add origin 5q and 8q resulting in duplications
7	Complex with markers	Concordant +	One mar is der(11)
8	-X and 6q deletion – 1 cell	Concordant +	t(12;21) and 14q deletion
9	Hyperdiploid, i(7q)	Concordant	None
10	Hyperdiploid, add(19p)	Concordant +	Add origin 12q
11	Hyperdiploid	Concordant +	+9 is der(9)t(8;9); +15 is +18
12	dic(9;20), +21	Concordant	None
13	t(1;19)(q23;p13.3)	Concordant +	14q deletion
14	t(1;19)(q23;p13.3)	Concordant	None
15	t(9;22)(q34;q11.2)	Concordant +	Homozygous 14q deletion
16	t(9;22)(q34;q11.2)	Concordant +	14q deletion
17	+i(21)(q10)x2	Concordant +	t(5;10), 16q and 22q deletions +14 is der(14)t(14;21), deletions (subtelomere) 18p and 15q
18	Hyperdiploid	Concordant +	t(12;21), 14q deletion
19	del(13)(q14q34)	Concordant +	del(3p) and 13 identified as der(3); deletion 16p and -20
20	Hypodiploid, t(9;22), dup(3q)	Missed duplication	t(12;21) and multiple rearrangements involving chromosomes 6, 9, 8, 10
21	Hyperdiploid with dup(1q)	Missed duplication	12p deletion
22	t(6;22)(q23;q12)	Concordant +	21q deletion, RUNX1 duplication
23	Complex with variant t(9;22)	Concordant +	21q deletion, 14q deletion, der(15)
24	add(21)(q22)	Concordant +	None
25	+X,+21c	Concordant +	Iso dic(17q)
26	del(12)(p11.2p13)	Concordant	None
27	Variant t(9;22)	Concordant +	None
28	del(6)(q13q27)	Concordant	Clarified add(10q) as iso dic(10q)
29	Complex with add(10q)	Concordant +	-19
30	Deletions 9p and 6q	Concordant +	Add origin 5q; der(19)t(1;19), 14q deletion
31	add(6q)	Concordant +	10p and 21q deletions
32	Hyperdiploid	Concordant +	None
33	9p (CDKN2A) deletion	Concordant	t(12;21), add material from chromosomes X and 11
34	12p deletion; add chromosome	Concordant+	t(12;21), 14q32 homozygous deletion, 8p deletion, 6q deletion, RUNX1 amplification
35	Complex; FISH positive for ETV6-RUNX1	Concordant+	CDKN2A dim (partial deletion), t(1;7)
36	t(7;17)	Concordant+	11q deletion; homozygous 14q deletion
37	t(2;14), +10	Concordant+	Homozygous CDKN2A deletion; 12p deletion, 14q deletion, 4p deletion
38	t(9;22)	Concordant+	Add material from chromosomes 13q, 3p, 11/9q
39	Complex karyotype; 'add' chromosomes	Concordant+	ICP 20 cells abnormal
40	Hyperdiploid; t(9;22); cyto 1 cell	Concordant+	mars are copies of extra chromosomes; TCF3 deletion (19p13.3) missed by ICP – not targeted
41	Hyperdiploid plus mars; FISH 19p deletion	Concordant+	6q deletion in 20 cells
42	6q deletion in two cells; 9p homozygous deletion	Concordant+	dic t(20;21) has tandem duplication of RUNX1 as well as telomere; 20q deletion confirmed
43	45,XY,dic(20;21)(q13.3;q22)	Concordant+	9p deletion; balanced t(5;11)
44	9q deletion in 80% of the cells	Concordant+	Marker chromosome, -12;
45	Marker chromosome, -12; ETV6-RUNX1 by FISH	Concordant+	t(12;19;21)
46	46,XX,t(2;7)(p21;p22)	Concordant+	14q deletion; possible t(2;12)
47	Hyperdiploid	Concordant+	Possible 12p deletion
48	Hyperdiploid	Concordant	None



### FIGURE LEGEND

A – C: Representative images showing a 3-fold increase in mitotic index; D: Trisomy 4 and 10; E: Normal chromosomes 2, 14 and 22 Mix 1; F: Normal chromosomes 8, 9 and 11 Mix 1; G: 14q32 (IGH) deletion; H: Homozygous 9p deletion and 6q deletion; I: der(6)(5;6)(q;q); J: RUNX1 Amplification; K – L: t(12;21)(p13;q22)

## DESIGN



## CHROMOSOME COLOR SCHEME

Color Pattern	Mix 1	Mix 2	Mix 3	Mix 4	Mix 5
Red	2	1	4	7	9p21.3
Green	14	12	10	16	6q21
R/G Hybrid	22	17	Y	13	13q14.2
Gold	8	X	5	3	21q22.12
Aqua	9	19	6	18	Xp22.33
Y/A Hybrid	11	21	20	15	

## CONCLUSION

- All numerical abnormalities identified by karyotype or FISH are detected by ALL-ICP.
- ALL-ICP refined the breakpoints and clarified marker and derivative chromosomes.
- Novel 14q32 (IGH) deletions and several previously unknown sub-telomere deletions were identified.
- There appears to be an association between t(9;22) and 14q deletions – nearly 50% with this deletion had the translocation.
- ALL common deletions and amplification of RUNX1 were detected.
- Non-recurring interstitial deletions and duplications not targeted by ALL-ICP were not detected.
- All known and new balanced translocations were detected.
- ALL-ICP detected additional abnormalities in 85% of the cases.
- A three-fold increase in mitotic index was observed in modified cultures.

ALL-ICP is a simple method that saves valuable time and detects all clinically relevant recurrent B- and T-cell chromosome abnormalities. Additionally, it characterizes the marker chromosomes and 'add' material. The improved diagnostic capabilities, supported by the study results, make ALL-ICP ideal for the genetic diagnostic workup of Acute Lymphocytic Leukemia on standard as well as modified cytogenetic harvests.

## REFERENCES & ACKNOWLEDGEMENTS

- Interphase Chromosome Profiling: A Method for Conventional Banded Chromosome Analysis Using Interphase Nuclei. Babu R, Van Dyke DL, Dev VG et al., Arch Pathol Lab Med. doi: 10.5858/ arpa.2016-0621-OA
- Clinical Significance of Cytogenetic Abnormalities in Adult Acute Lymphoblastic Leukemia. Faderl S, Kantarjian HM, Talpaz M, and Estrov Z. Blood, Vol 91, No 11 (June 1), 1998: pp 3995-4019
- Cytogenetics and Molecular Genetics of Acute Lymphoblastic Leukemia. Mrózek K, Harper DP, and Aplan PD. Hematol Oncol Clin North Am. 2009 October ; 23(5): 991–v. doi:10.1016/j.hoc.2009.07.001
- Microdeletions are a general feature of adult and adolescent acute lymphoblastic leukemia: Unexpected similarities with pediatric disease. Paulsson K, Cazier JB, Finlay MacDougall F et al., PNAS, vol. 105, no. 18 , May 6, 2008: 6708–6713

We would like to thank Yvonne Banol (InteGen LLC) for technical support.

Please contact us at [contact@integenllc.com](mailto:contact@integenllc.com)  
[www.integenllc.com](http://www.integenllc.com)

